

Product Datasheet

TMSB15A antibody (orb18498)

Description	Goat polyclonal antibody to TMSB15A
Species/Host	Goat
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA
Target	Thymosin-like 8
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:2000.WB: No signal obtained yet but low background observed in 293 and A549 lysates at up to 2 µg/ml (please note that our SDS PAGE/western blotting protocol is not optimized for detection of proteins smaller than 15kDa). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates
Protein Sequence	SDKPDLSEVEKFDR
Clonality	Polyclonal
MW	5.23
Entrez	11013
Dilution Range	ELISA: 1:2000, WB: 2 µg/ml
Expiration Date	12 months from date of receipt.