

Product Datasheet

TRIM35 antibody (orb18499)

Description

Goat polyclonal antibody to TRIM35

Species/Host

Goat

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, WB

Applications
Target

TRIM35 / MAIR

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

ELISA: Peptide ELISA: antibody detection limit dilution 1:32000. WB: Approx. 25-30kDa band observed in human testis and human brain lysates at 1ug of primary, this band was successfully blocked by incubation with the immunising peptide. Please note that the band we observe is different from the predicted sizes of approx. 57kDa and 24kDa according to NP_055881 and NP_741983. Currently we cannot find an explanation for this in the literature. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates Have any splice variants/modified forms been reported

Protein Sequence

ERSPDVSPGPSRS

Clonality

Polyclonal

MW


56.5

Entrez
23087
Dilution Range

ELISA: 1:32000, WB: 1-3 µg/ml

Expiration Date

12 months from date of receipt.



Western blot analysis of K562 lysate using...