



Product Datasheet

RPS27 antibody (orb18500)


Description	Goat polyclonal antibody to RPS27
Species/Host	Goat
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	MPS1 / RPS27
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:64000.WB: Approx 12kDa band observed in lysates of T-cell line Jurkat (calculated MW of 9.5kDa according to NP_001021.1). Recommended concentration: 0.05-0.2 µg/ml.
Protein Sequence	PLAKDLLHPSPEE
Clonality	Polyclonal
MW	9.5
Entrez	6232
Dilution Range	ELISA: 1:128000, WB: 0.03-0.1 µg/ml, IHC-P: 3.75 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Jurkat lysate u...



Western blot analysis of Jurkat lysate u...



Immunohistochemical staining of Human To...