



Product Datasheet

LSM12 antibody (orb185079)

Description	LSM12 antibody
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human LSM12 (101-195/195 aa)
Target	LSM12
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	22 kDa
Uniprot ID	Q3MHD2
Dilution Range	IHC-P=1:100-500
Expiration Date	12 months from date of receipt.



IHC-P
analysis of
mouse testis
using LSM...



IHC-P
analysis of
rat liver using
LSM12 ...