


Product Datasheet

Rat ROR1 Protein (orb18509)

Description	Tyrosine-protein kinase transmembrane receptor ROR1 is also known as Neurotrophic tyrosine kinase, r...
Reactivity	Rat
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	ROR1
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 43.5 kDa. The protein migrates as 60-67 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	EDL97819
Purity	95%
Source	Rat ROR1, His Tag (orb18509) is expressed from human 293 cells (HEK293). It contains AA Gln 30 - Glu 403 (Accession # EDL97819).
MW	43.5 kDa
NCBI	EDL97819
Expiration Date	6 months from date of receipt.



SDS-PAGE
analysis of
Rat ROR1
protein