



Product Datasheet

Human VEGF R3 / FLT4 Protein (orb18510)

Description	Vascular endothelial growth factor receptor 3 (VEGF R3), also known as FLT-4, together with the othe...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	VEGF R3
Preservatives	Tris with Glycine, Arginine and NaCl, pH7.5
Form/Appearance	Powder
Storage	-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.
Tag	C-Fc
Note	For research use only
Application notes	This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 111.2 kDa. As a result of glycosylation, the protein migrates as 140-160 kDa under reducing (R) condition, and 280-320 kDa under non-reducing (NR) condition (SDS-PAGE).
Purity	95%
Source	Human VEGF R3, Fc Tag (orb18510) is expressed from human 293 cells (HEK293). It contains AA Tyr 25 - Ile 776 (Accession # NP_002011.2).
MW	111.2 kDa
NCBI	NP_002011.2
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Human VEGF R3 prote...



Human VEGF R3, Fc Tag (Cat. No. orb18510...