

Product Datasheet

TRAIP antibody (orb18513)

Description

Goat polyclonal antibody to TRAIP


Species/Host

Goat

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, IHC

Applications

Immunohistochemical
staining of Human
Sp...

Target

TRIP / TRAIP

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

ELISA: Peptide ELISA: antibody detection limit dilution 1:2000.WB: Approx. 40kDa band observed in Human Testis lysate at 1ug of primary. Please note that the band we observe is different from the predicted size of approx. 54kDa according to NP_005870. Currently we cannot find an explanation for this in the literature. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates Have any splice variants/modified forms been reported

Protein Sequence

SLFQAKLDTFLWS

Clonality

Polyclonal

MW

53.3

Entrez

10293

Dilution Range

ELISA: 1:1000, IHC-P: 4-6µg/ml

Expiration Date

12 months from date of receipt.