



Product Datasheet

GPS1 antibody (orb18514)


Description	Goat polyclonal to GPS1
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	GPS1 / COPS1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	This antibody is expected to recognise both reported isoforms.ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 55kDa band observed in Human Testis lysate (predicted MW of 58kDa according to NP_004118). Recommended for use at 0.5 µg/ml.
Protein Sequence	PANSQSRMSTNM
Clonality	Polyclonal
MW	59.1; 55.5
Entrez	2873
Dilution Range	ELISA: 1:32000, WB: 0.5 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



orb18514 staining
(0.5ug/ml) of Human
Te...



Immunohistochemical
staining of Human
Pr...



Immunohistochemical
staining of Human
Te...