


Product Datasheet

TRIM4 antibody (orb18515)

Description	Goat polyclonal to TRIM4
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	TRIM4
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	This antibody is expected to recognise both reported human isoforms (as represented by NP_148977 and NP_149082).ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 55-60kDa band observed in HepG2 whole cell lysate (predicted MW of 58kDa and 55kDa according to NP_148977 and NP_149082). Recommended for use at 1 µg/ml.
Protein Sequence	ASLVIPPVTDRK
Clonality	Polyclonal
MW	57.5, 54.5
Entrez	89122
Dilution Range	ELISA: 1:32000, WB: 1 µg/ml
Expiration Date	12 months from date of receipt.



orb18515
staining
(1ug/ml) of
HepG2 lysa...