



Product Datasheet

Myosin-3 antibody (orb185150)

Description	Myosin-3 antibody
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human Myosin-3 (1801-1940/1940 aa)
Target	Myosin-3
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	224 kDa
Uniprot ID	P11055
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Rat skeletal muscle tissue.



Immunohistochemical staining of Mouse placenta tissue.