


Product Datasheet

USP4 antibody (orb18517)

Description	Goat polyclonal antibody to USP4
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	USP4 / UNP
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:8000.WB: No signal obtained yet but low background observed in Human Brain lysate at up to 1 µg/ml. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates
Protein Sequence	AEGGGCRERPDAE
Clonality	Polyclonal
MW	109; 104
Entrez	7375
Dilution Range	ELISA: 1:8000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Heart lys...