

Product Datasheet

MAPK8IP3 antibody (orb18518)

Description	Goat polyclonal antibody to MAPK8IP3
Species/Host	Goat
Reactivity	Bovine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Target	JIP3 / Syd2 / JSAP1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:4000.WB: No signal obtained yet but low background observed in Human Heart lysate at up to 2 µg/ml. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates
Protein Sequence	MEIQMDEGGGVV
Clonality	Polyclonal
MW	148, 147
Entrez	23162
Dilution Range	ELISA: 1:4000, WB: 2 µg/ml, IHC-P: 4-6µg/ml
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Mouse Br...