



Product Datasheet

MLH3 antibody (orb18532)

Description	Goat polyclonal antibody to MLH3
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF
Target	MLH3
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:2000.WB: No signal obtained yet but low background observed in HeLa lysate at up to 2 µg/ml. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates
Protein Sequence	RQSLQQSMPPCEPP
Clonality	Polyclonal
MW	164; 161
Entrez	27030
Dilution Range	ELISA: 1:2000, IHC-P: 10ug/ml
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of paraforma...



Flow cytometric analysis of paraformalde...