


Product Datasheet

USP15 antibody (orb18533)

Description	Goat polyclonal antibody to USP15
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Target	USP15
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:16000.WB: No signal obtained yet but low background observed in Human Skeletal Muscle and Kidney extracts at up to 1 µg/ml. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates
Protein Sequence	NDIENENCMHTN
Clonality	Polyclonal
MW	109
Entrez	9958
Dilution Range	ELISA: 1:2000, IHC-P: 2.5-3.8µg/ml
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Human Liver tissue