

Product Datasheet

MAPK15 (phospho-T175+Y177) antibody (orb185465)

Description	MAPK15 (phospho-T175+Y177) antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthesised phosphopeptide derived from human MAPK15 around the phosphorylation site of Thr175 + Tyr177 (AV(p-T)E(p-Y)VA)
Target	phospho-MAPK15 (Thr175 + Tyr177)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	60 kDa
Uniprot ID	Q8TD08
Expiration Date	12 months from date of receipt.