



Product Datasheet

MCL1 (phospho-Ser159) antibody (orb185519)

Description	MCL1 (phospho-Ser159) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ICC, IHC-P, WB
Immunogen	KLH conjugated synthesised phosphopeptide derived from human Mcl 1 around the phosphorylation site of Ser159. (DG(p-S)LP)
Target	phospho-MCL1 (Ser159)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	37 kDa
Uniprot ID	Q07820
Dilution Range	WB=1:500-2000, IHC-P=1:100-500, ICC=1:100
Expiration Date	12 months from date of receipt.



ICC analysis of A549 cell using phospho-...



IHC-P analysis of mouse brain using phos...