



Product Datasheet

RANBP9 antibody (orb18558)

Description	Goat polyclonal antibody to RANBP9
Species/Host	Goat
Reactivity	Canine, Human, Mouse, Porcine
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IHC
Target	RANBPM / RANBP9
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:16000.IHC: In paraffin embedded human tonsil, stains nuclei in epithelial cells of tonsillar cryptae. Recommended concentration, 3ug/ml.WB: Preliminary experiments gave no signal but low background in whole cell Hela lysates at up to 1 µg/ml.
Protein Sequence	GSCAFATVEDYLH
Clonality	Polyclonal
MW	77.8; 55.1
Entrez	10048
Dilution Range	ELISA: 1:16000, WB: 1 µg/ml, IHC-P: 10 µg/ml
Expiration Date	12 months from date of receipt.



orb18558 (3ug/ml)
staining of paraffin
e...



Immunohistochemical
staining of Human
Te...