



Product Datasheet

SIL1 antibody (orb18559)


Description	Goat polyclonal antibody to SIL1
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	BAP / SIL1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:8000.WB: Approx 55-60 kDa band observed in Human Liver and Kidney extracts (predicted MW of 53kDa according to NP_071909.1). Recommended for use at 1-2 µg/ml.
Protein Sequence	ELLGSVNSLLKELR
Clonality	Polyclonal
MW	52.6
Entrez	64374
Dilution Range	ELISA: 1:8000, WB: 1-2 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Kidney ly...



Immunohistochemical staining of Human Te...



Immunohistochemical staining of Human li...