


Product Datasheet

ZNF804A antibody (orb1856)

Description	ZNF804A antibody
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human ZNF804A (75-170/1209 aa)
Target	ZNF804A
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	137 kDa
Uniprot ID	Q7Z570
Dilution Range	IHC-P=1:100-500
Expiration Date	12 months from date of receipt.



Paraffin-
embedded
rat brain
tissue fixed...