

Product Datasheet

OXSR1 (phospho-Thr185) antibody (orb185728)

Description	OXSR1 (phospho-Thr185) antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	KLH conjugated synthesised phosphopeptide derived from human OXSR1 around the phosphorylation site of Thr185 (RK(p-T)FV)
Target	phospho-OXSR1 (Thr185)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	58 kDa
Uniprot ID	O95747
Expiration Date	12 months from date of receipt.