

Product Datasheet

MTCH1 antibody (orb18577)

Description	Goat polyclonal antibody to MTCH1
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA
Target	MTCH1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:16000.WB: Preliminary experiments gave no signal but low background in Human Heart and Placenta lysates at up to 1.5 µg/ml. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates
Protein Sequence	FRRVSSGSCFALE
Clonality	Polyclonal
MW	39.9
Entrez	23787
Dilution Range	ELISA: 1:16000, WB: 1.5 µg/ml
Expiration Date	12 months from date of receipt.