

Product Datasheet

SHFM1 antibody (orb18579)

Description

Goat polyclonal antibody to SHFM1

Species/Host

Goat

Reactivity

Bovine, Canine, Human, Mouse, Porcine

Conjugation

Unconjugated

Tested

ELISA, ICC

Applications
Target

DSS1 / SHFM1

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

ELISA: Peptide ELISA: antibody detection limit dilution 1:32000. WB: No signal obtained yet but a preliminary experiment showed low background in Human Liver extract at up to 1 µg/ml. This used our standard western blotting protocol which we would not expect to detect proteins as small as the predicted size of approx. 8KDa according to NP_006295. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates

Protein Sequence

RAELEKHGYKMETS

Clonality

Polyclonal

MW

8.3

Entrez

7979

Dilution Range

ELISA: 1:32000, WB: 1 µg/ml

Expiration Date

12 months from date of receipt.