

Product Datasheet

ORMDL1 + ORMDL2 + ORMDL3 antibody [Out of stock] (orb185868)

Description	ORMDL1 + ORMDL2 + ORMDL3 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Immunogen	KLH conjugated synthetic peptide derived from human ORMDL1 + ORMDL2 + ORMDL3 (1-80/153 aa)
Target	ORMDL1 + ORMDL2 + ORMDL3
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	17 kDa
Uniprot ID	Q53FV1, Q9P0S3, Q8N138
Expiration Date	12 months from date of receipt.