

Product Datasheet

phosphor-PPP3CA (phospho-Ser469) antibody (orb185888)

Description	phosphor-PPP3CA (phospho-Ser469) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ICC, IHC-P
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human PPP3CA around the phosphorylation site of Ser469 (421-521/521 aa)
Target	Phosphor-PPP3CA (Ser469)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	57 kDa
Uniprot ID	Q08209
Dilution Range	IHC-P=1:100-500, ICC=1:25
Expiration Date	12 months from date of receipt.