



Product Datasheet

BAG3 antibody (orb18594)


Description	Goat polyclonal antibody to BAG3
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	BAG3 / BIS/ CAIR1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:8000.WB: Approx 85kDa band observed in Human Skeletal Muscle lysates (calculated MW of 61.6kDa according to NP_004272.2). The observed molecular weight corresponds to earlier findings in literature with different antibodies (Iwasaki et al, Cancer Res. 2007 Nov 1;67(21):10252-9. PMID: 17974966). Recommended concentration: 0.1-0.3 µg/ml.
Protein Sequence	SSMTDTPGNPAAP
Clonality	Polyclonal
MW	61.6
Entrez	9531
Dilution Range	ELISA: 1:2000, WB: 0.3-1 µg/ml, IHC-P: 3.75 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Skeletal ...



Immunohistochemical staining of Human He...



Immunohistochemical staining of Human Sk...