



Product Datasheet

pIRF3 antibody (orb1860)


Description	pIRF3 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	KLH conjugated synthesised phosphopeptide derived from human IRF3 around the phosphorylation site of Ser386 (S(p-S)LE)
Target	phospho-IRF3 (Ser386)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	47 kDa
Uniprot ID	Q14653
Dilution Range	WB=1:500-2000, IHC-P=1:100-500, Flow-Cyt=1ug/test, IHC-F=1:100-500, ICC=1:100, IF=1:50-200, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Human Liver



Immunohistochemical staining of Human Gastric



Western blot analysis of mouse bone lysate