



Product Datasheet

SOCS6 antibody (orb186020)

Description	SOCS6 antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human SOCS6 (21-120/535 aa)
Target	SOCS6
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	59 kDa
Uniprot ID	O14544
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



Western blot:Heart
(Mouse);Primary:Anti-...



Western blot:Heart
(Chicken);Primary:Ant...