


Product Datasheet

SP1 (phospho-Thr739) antibody (orb186031)

Description	SP1 (phospho-Thr739) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC
Immunogen	KLH conjugated synthesised phosphopeptide derived from human SP1 around the phosphorylation site of Thr739 (TA(p-T)PS)
Target	phospho-SP1 (Thr739)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	81 kDa
Uniprot ID	P08047
Dilution Range	ICC=1:100
Expiration Date	12 months from date of receipt.



ICC analysis of A549 cell using phospho-...