


Product Datasheet

XPO7 antibody (orb18605)

Description	Goat polyclonal antibody to XPO7
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	RANBP16 / Exportin 7
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 110kDa band observed in Human Testis and Hela lysates (predicted size of approx. 124kDa according to NP_055839). Recommended for use at 0.5-1 µg/ml
Protein Sequence	NSTYGVNSNDMMMS
Clonality	Polyclonal
MW	125; 124; 124
Entrez	23039
Dilution Range	ELISA: 1:32000, WB: 0.5-1 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Human Testis tissue



Western blot analysis of Hela lysate using XPO7 antibody