


Product Datasheet

Scavenger Receptor BI + BII antibody (orb186078)

Description	Scavenger Receptor BI + BII antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	KLH conjugated synthetic peptide derived from human Scavenger Receptor BI + BII (31-130/552 aa)
Target	Scavenger Receptor BI + BII
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	61 kDa
Uniprot ID	Q8WTV0
Dilution Range	WB=1:500-2000, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse adrenal g...