



Product Datasheet

KPNA4 antibody (orb18608)


Description	Goat polyclonal antibody to KPNA4
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	KPNA4 / IPOA3
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.IHC: In paraffin embedded Human Testis, shows strong nuclear staining of some cells in the seminiferous tubules. Recommended concentration, 3ug/ml. In paraffin embedded Human Small Intestine shows pixulate staining in nuclei of enterocytes. Recommended concenWB: Approx 55-60kDa band observed in Human Testis and Ovary lysates (predicted size of approx. 58kDa according to NP_002259.1). Recommended for use at 0.2-1 µg/ml
Protein Sequence	NSSANVPTEGFQF
Clonality	Polyclonal
MW	57.9
Entrez	3840
Dilution Range	ELISA: 1:128000, WB: 0.5-1 µg/ml
Expiration Date	12 months from date of receipt.



orb18608 staining (0.5ug/ml) of Human Te...



orb18608 (3ug/ml) staining of paraffin e...



orb18608 (5ug/ml) staining of paraffin e...