



Product Datasheet

KPNA3 antibody (orb18610)

Description	Goat polyclonal antibody to KPNA3
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, IHC, WB
Target	KPNA3 / IPOA4
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.IHC: In paraffin embedded Human Heart shows strong pixulate staining of myocardial fibres in transverse section. Recommended concentration, 5-10ug/ml.WB: Approx 55-60kDa band observed in Human Testis and Kidney lysates (predicted size of approx. 58kDa according to NP_002258). Recommended for use at 0.1-0.5 µg/ml
Protein Sequence	DPTANLQTKEFNF
Clonality	Polyclonal
MW	57.8
Entrez	3839
Dilution Range	ELISA: 1:128000, WB: 0.1-0.3 µg/ml
Expiration Date	12 months from date of receipt.



orb18610
staining
(0.1ug/ml) of
Human Te...



orb18610
(5ug/ml)
staining of
paraffin e...



Western blot
analysis of
HeLa lysate
usi...