

Product Datasheet

NUPR1 antibody (orb18613)

Description

Goat polyclonal antibody to NUPR1

Species/Host

Goat

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA

Applications
Target

P8 / NUPR1

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

ELISA: Peptide ELISA: antibody detection limit dilution 1:32000. WB: Preliminary experiments gave no signal but low background in HepG2 and Human Pancreas extracts at up to 1 µg/ml. Please note that we would not expect our routine western blot protocol to reliably detect proteins as small as the predicted size of this protein (approx 9kDa according to NP_036517). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates

Protein Sequence

ATFPATSAPQQP

Clonality

Polyclonal

MW

10.7; 8.9

Entrez

26471

Dilution Range

ELISA: 1:32000, WB: 1 µg/ml

Expiration Date

12 months from date of receipt.