

Product Datasheet

Simian Rotavirus VP4 antibody (orb186144)

Description	Simian Rotavirus VP4 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from Simian Rotavirus VP4 (201-300/776 aa)
Target	Rotavirus A VP4
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	58/85 kDa
Uniprot ID	P12473
Expiration Date	12 months from date of receipt.