

Product Datasheet

ATP6AP2 antibody (orb18618)

Description	Goat polyclonal antibody to ATP6AP2
Species/Host	Goat
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, WB
Target	ATP6IP2 / Renin receptor
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:16000.WB: Approx 37kDa band observed in Human Kidney lysates (calculated MW of 39.0kDa according to NP_005756.2). Recommended for use at 0.25-1 µg/ml
Protein Sequence	SIYRMTNQKIRMD
Clonality	Polyclonal
MW	39
Entrez	10159
Dilution Range	ELISA: 1:16000, WB: 1-2 µg/ml, IHC-P: 10ug/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Kidney ly...