

Product Datasheet

AEV polyprotein antibody (orb1862)

Description	AEV polyprotein antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	KLH conjugated synthetic peptide derived from AEV polyprotein, RNA-directed RNA polymerase 3D-POL (2001-2134/2134 aa)
Target	AEV polyprotein
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9YLS4
Expiration Date	12 months from date of receipt.