



Product Datasheet

SRGAP1 antibody (orb186233)

Description	SRGAP1 antibody
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human SRGAP1 (51-150/1085 aa)
Target	SRGAP1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	124 kDa
Uniprot ID	Q7Z6B7
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of rat testis...



Immunohistochemical staining of rat brain...