


Product Datasheet

EIF4ENIF1 antibody (orb18625)

Description	Goat polyclonal antibody to EIF4ENIF1
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	4E-T / EIF4ENIF1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.IHC: In paraffin embedded Human Pancreas shows strong staining of cytoplasm in acinar cells-. Recommended concentration, 5-10ug/ml.WB: Approx 140-150kDa band observed in 293 lysates (predicted size of approx. 108kDa according to NP_062817.1, however our observation agrees with that of Dostie et al (see below)). Recommended for use at 0.25-0.5 µg/ml
Protein Sequence	AKVISVDELEYRQ
Clonality	Polyclonal
MW	108; 88.2
Entrez	56478
Dilution Range	ELISA: 1:32000, WB: 0.25-0.5 µg/ml, IHC-P: 5-10µg/ml
Expiration Date	12 months from date of receipt.



orb18625
staining
(0.25ug/ml)
of 293 lys...



orb18625
(5ug/ml)
staining of
paraffin e...