

Product Datasheet

Troponin T-Slow Skeletal antibody (orb186263)

Description	Troponin T-Slow Skeletal antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human Troponin T-Slow Skeletal (61-160/278 aa)
Target	Troponin T-Slow Skeletal
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	33 kDa
Uniprot ID	P13805
Expiration Date	12 months from date of receipt.