

Product Datasheet

Troponin I fast + slow skeletal muscle antibody (orb186265)

Description	Troponin I fast + slow skeletal muscle antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human TNNI1+TNNI2 (51-150/182 aa)
Target	Troponin I fast + slow skeletal muscle
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	33 kDa
Uniprot ID	P19237
Expiration Date	12 months from date of receipt.