


Product Datasheet

Troponin T fast + slow skeletal muscle antibody (orb186267)

Description	Troponin T fast + slow skeletal muscle antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human Troponin T fast + slow skeletal muscle (31-130/269 aa)
Target	Troponin T fast + slow skeletal muscle
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	32 kDa
Uniprot ID	P45378
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse skin lysa...