



Product Datasheet

TAGLN antibody (orb18630)


Description	Goat polyclonal to TAGLN
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	SM22 Alpha / Transgelin
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	NP_003177.2 and NP_001001522.1 represent identical protein.ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 24kDa band observed in Human Placenta lysates (calculated MW of 22.6kDa according to NP_003177.2 and NP_001001522.1). Recommended concentration: 0.01-0.03 µg/ml.
Protein Sequence	MTGYGRPRQIIS
Clonality	Polyclonal
MW	22.6
Entrez	6876
Dilution Range	ELISA: 1:4000, WB: 0.1-1 µg/ml, IHC-P: 2.5µg/ml
Expiration Date	12 months from date of receipt.



orb18630 (0.03µg/ml)
staining of Mouse ...



orb18630 (0.03µg/ml)
staining of Human ...



Immunohistochemical
staining of Human
Pr...