

Product Datasheet

TRIM29 (acetyl K116) antibody (orb186338)

Description	TRIM29 (acetyl K116) antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated Synthesised acetylpeptide derived from human TRIM29 around the acetylation site of K116 (AK(Ac-K)PP)
Target	TRIM29 (acetyl K116)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	66 kDa
Uniprot ID	Q14134
Expiration Date	12 months from date of receipt.