



Product Datasheet

RANGAP1 antibody (orb18635)

Description	Goat polyclonal antibody to RANGAP1
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	RANGAP1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.IHC: In paraffin embedded Human Liver shows vesicular staining in the cytoplasm of hepatocytes. Recommended concentration, 2-4ug/ml.WB: Approx 70kDa band observed in Human Brain and 3T3 lysates (predicted MW of 68kDa according to NP_002874). Recommended for use at 0.1-0.5 µg/ml.
Protein Sequence	ASEDIAKLAETLAK
Clonality	Polyclonal
MW	63.5
Entrez	5905
Dilution Range	ELISA: 1:32000, WB: 0.1-0.5 µg/ml, IHC-P: 2-4 µg/ml
Expiration Date	12 months from date of receipt.



orb18635
staining
(0.2ug/ml) of
3T3 lysa...



orb18635
(2.5ug/ml)
staining of
paraffin...