

Product Datasheet

SNTG1 antibody (orb18638)

Description	Goat polyclonal antibody to SNTG1
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Target	Syntrophin gamma 1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:8000.WB: Preliminary experiments gave no signal but low background in Human Testis and Brain extracts at up to 1 µg/ml. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates
Protein Sequence	ASSATTSKAKYTT
Clonality	Polyclonal
MW	58
Entrez	54212
Dilution Range	ELISA: 1:8000, WB: 1 µg/ml
Expiration Date	12 months from date of receipt.