


Product Datasheet

HTATIP2 antibody (orb18640)

Description	Goat polyclonal antibody to HTATIP2
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	TIP30 / CC3
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 28kDa band observed in lysates of cell line HeLa (calculated MW of 27.0kDa according to NP_006401.3). Recommended concentration: 0.3-1 µg/ml.
Protein Sequence	AETEALSKLREDFR
Clonality	Polyclonal
MW	30.1; 27.0; 14.9
Entrez	10553
Dilution Range	ELISA: 1:32000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.



orb18640
(0.6ug/ml)
staining of
HeLa lys...