

Product Datasheet


MYRIP antibody (orb18642)

Description	Goat polyclonal antibody to MYRIP
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	MYRIP / SLAC2C
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:64000.WB: Approx 90kDa band observed in Human Brain (Amygdala) lysates (calculated MW of 95.7kDa according to NP_056275.2). Recommended concentration: 0.3-1 µg/ml.
Protein Sequence	KDLMEPALESVMY
Clonality	Polyclonal
MW	95.7
Entrez	25924
Dilution Range	ELISA: 1:64000, WB: 0.3-1 µg/ml, IHC-P: 3.75 µg/ml
Expiration Date	12 months from date of receipt.




Western blot analysis of Human Amygdala

...



Immunohistochemical staining of Human Br...



Immunohistochemical staining of Human Co...