



Product Datasheet

MAP4K1 antibody (orb18655)

Description	Goat polyclonal antibody to MAP4K1
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	HPK1 / MAP4K1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 95-100kDa band observed in Jurkat (predicted MW of 96kDa according to NP_009112). Recommended for use at 0.5-2 µg/ml.
Protein Sequence	DVVDPDFIFNRDPR
Clonality	Polyclonal
MW	90.4; 91.3
Entrez	11184
Dilution Range	ELISA: 1:32000, WB: 0.5-2 µg/ml, IHC-P: 2.5 µg/ml
Expiration Date	12 months from date of receipt.



orb18655 staining
(0.5ug/ml) of Jurkat
l...



Immunohistochemical
staining of Human
Th...