



Product Datasheet

PPP2R5A antibody (orb18661)

Description	Goat polyclonal antibody to PPP2R5A
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	PPP2R5A
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 55kDa band observed in human heart and skeletal muscle lysate (calculated MW of 56.2kDa according to NP_006234.1). Recommended for use at 0.5-2 µg/ml.
Protein Sequence	AYNMHSILSNTSAE
Clonality	Polyclonal
MW	56.2
Entrez	5525
Dilution Range	ELISA: 1:32000, WB: 0.5-2 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of human muscle ly...



Immunohistochemical staining of Human He...