


Product Datasheet

VDAC2 antibody (orb18662)


Description	Goat polyclonal antibody to VDAC2
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	VDAC2 (C Terminus)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Concentration	1-3 µg/ml
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:64000.
Protein Sequence	GHKVGLALELEA
Clonality	Polyclonal
MW	31.6; 33.3; 27.4
Entrez	7417
Dilution Range	ELISA: 1:16000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.



orb18662
Negative
Control
showing
staini...



orb18662
(6µg/ml)
staining of
paraffin ...



orb18662
(1ug/ml)
staining of
CaCo-2 (A)...